

→ 1. Set patient up


Advise the patient to download the SleepCare App and register.

Tip – search 'SleepCare App' to find it
or use this link – [SleepCare – Good Sleep Co](#)


The rest of the patient set up (B2C) is outlined here.



Apple



Android



SleepCare
APP



Getting Started (Health Provider's Viewpoint)

→ 2. Select Test Type

- Good Sleep Test, or
- Good Sleep Test Specialist Reported

→ 3. Patient Pricing


Patients can acquire a 1 Good Sleep Test device that can be used for 90 days worth \$50.

A specialist-reported test provides a report and diagnoses from a Respiratory and Sleep Physician and costs clinicians an additional \$50 per report. This option will also allow for private health rebates should the patient require treatment for sleep disordered breathing.



Monitoring using the Good Sleep Test

→ 4. Create a Test



→ 5. Start a Test

Click Drop Down Arrow and Choose 'Good Sleep Test' or 'Good Sleep Test Specialist Reported'.

Note for Specialist Reported you need to answer General and STOP BANG questions and +\$50. Indicate nights and Gender.



If you would like more information,
please contact us:

hello@goodsleepco.health



Monitoring using the Good Sleep Test

→ 6. Test results

In the Test board, you will see the status if it is 'Pending' (Test is not yet done) and 'Completed'.

Results will be seen beside the 'Completed' Status.

DATE ORDERED	SLEEP STUDY	PATIENT	NIGHTS ORDERED	STATUS
30 May 2024	Good Sleep Test	Bill Johnson	1	Pending
28 May 2024	Good Sleep Test	Jane Smith	1	Completed

Troubleshooting Tips

→ Patient birth date correct

Ensure date of birth of patient is correct when creating the patient – this cannot be updated afterwards.

→ Control Bluetooth

Ensure there are no other Bluetooth devices in the room while doing the test

→ Check connection

Verify that you have a good internet connection by navigating to any webpage in your browser, for example:

- If the Good Sleep Test device isn't visible, try refreshing the app by logging out and then logging back in.
- If you are on Wifi, try switching to cellular (or vice-versa)
- If you have a VPN configured on your phone, that may cause issues. If you can disable it, try that.
- If you have another mobile device in the household, such as a family members' phone or a tablet, try that and see if you are successful there.
- If the recording time or quality is low once test results are received, please contact us hello@goodsleepco.health



If you would like more information,
please contact us:

hello@goodsleepco.health