



Good Sleep Test

GoodSleep^{co}

Overarching brand that describes
our purpose

hushd[™]

Our oral appliances

The logo for Resdent features a stylized blue flame icon above the word "Resdent" in a blue sans-serif font.

The system that manages the tools needed for HCP's
to operate in the dental sleep medicine space.

The logo for SleepCare features a stylized blue flame icon above the word "SleepCare" in a blue sans-serif font.

POWERED BY RESDENT

The app that takes care of our
patients. Like a practitioner in their
pocket

What is the Good Sleep Test

- 1. Level 4 testing – Mechanism**
- 2. The Good Sleep Test hardware**
- 3. How to set up a test and ResDent**
- 4. The Good Sleep Test in the SleepCare app**
- 5. The Good Sleep Test reports**



1. Level 4 testing

Brand of level 4 testing

- WatchPat
- NightOwl
- Sunrise
- Sleep Image
- **Good Sleep Test**

How do they work

Pulse Oximetry

PPG

HR

Actigraphy

Pulse Ox

- Helps us understand desaturation of oxygen from apnea, hypopnea and RERA (Respiratory effort related arousal – meaning you're trying to breathe and its not going to plan and your brain creates an arousal in response)
- Average SPO2 and CV risk
- Highs and lows help us assess risk

PPG

A reliable method for determining sympathetic activation – AKA sleep arousals or arousals associated with apnea, hypopnea or RERA

We also use the PPG to help with determining sleep stage in combination with actigraphy

PPG can measure arterial contraction with respiratory sleep disturbances - Metaboreflex

Heart rate

With sympathetic arousal we usually see HR increase.

We use this to determine whether an SB event has occurred

Actigraphy

- Are you moving in your sleep and when
- We use this to determine the stage of sleep
- In combination with the other variables this gives us the last piece of information we need to determine the stage of sleep
- We need to know deep non REM, REM and light non REM to determine how OSA, UARS (usually diagnosed from RDI aka RERA) and how much the sleep breathing problem is affecting sleep architecture
- Poor sleep architecture = patient feels like a sack of potatoes



2. The Good Sleep Test Hardware and product



good
sleep^{TEST}

good sleep^{TEST}
SPECIALIST REPORT

Patient
Name: Phil (45) male
DOB: 5 Jan 1978 (26 years)
History: Patient is a 26-year-old male with a BMI of 26.5. A routine sleep study was conducted to evaluate for sleep apnoea.

Test
Device: Sleepio
Date: 6 Jun 2023

Physician
Name: Dr. David Smith
GoodSleepco (Pty) Ltd
Certified Sleep Medicine for the SBA
Provider No: 2024/001

Interpretation
On the night of 2023/06/04, patient completed an unattended home sleep test (HST) evaluation. During the night, 100% of channels were recorded and all data was available for analysis. The following table provides a summary of the results. The following is the interpretation for this study night.

Impression
Moderate to severe Sleep Apnoea/Respiratory Disturbance

Recommendation
Further assessment and treatment of Sleep Apnoea/Respiratory Disturbance is recommended. A referral to a Sleep Medicine Specialist is recommended. Further assessment and treatment of Sleep Apnoea/Respiratory Disturbance is recommended. A referral to a Sleep Medicine Specialist is recommended.

Results
AAS: 26.5
AHI: 29.6
94% SPO2

GoodSleepco
Patient Name: Phil - 1 Data 2023



Good Sleep Test

GoodSleepco

good
sleep^{TEST}

HOME SLEEP TEST



TGA
Approved

Multi-night Testing
Obstructive & Central Apnea Diagnosis
Specialist Reported Sleep Studies
(TO HELP WITH PRIVATE INSURANCE REIMBURSEMENT
OF CPAP AND ORAL APPLIANCE THERAPY)

Simple. Accurate. Cost effective.

GoodSleep^{co}

goodsleepco.health



3. How to set up a test and ResDent

Setting up testing

- Autoscored studies
- Multiple night studies
- Adding a sleep physician test
- Reports generated for practitioner/pharmacist/pharmacy assistance and possibly patient if selected



5. The Good Sleep Test in the SleepCare app

What it looks like for the patient

.

- <https://xd.adobe.com/view/9cc11583-c36d-4cc3-877d-ddb4b20aed00-c992/screen/9d397bae-e79b-4e2d-9383-791dfdbcb0ce/>



6. The Good Sleep Test reports

Patient

Name: Paul McKenzie

DOB: 3 Jan 1987 (36 years)

History: Patient is a 36 year old Male with a BMI of 26.5. A home sleep study was performed to evaluate for sleep apnea.

Test

Device: Aurora

Date: 6 Jun 2023

Physician

Name: Howard Dean, MD

Qualification: Board Certified In Sleep Medicine by the ABIM

Provider No. 269410DJ

Interpretation

On the night of 2019/10/14, patient completed an unattended home sleep test (HST) evaluation (Aurora, type IV HST), with channels monitoring pulse, oxygen saturation and respiratory analysis by peripheral arterial tonometry. Total sleep time was determined through analysis of photoplethysmography signals. Sleep time was determined by classifying each epoch as "sleep" or "wake" and summing epochs of sleep. The following is the interpretation for that study night.

Impression

Moderate Obstructive Sleep Apnea/Hypopnea Syndrome.

Recommendation

Treatment options include oral appliance or CPAP therapy further pending patient preference. If amenable to PAP therapy, recommend APAP 5-15 cmH2O.

Given severity, recommend CPAP the gold standard therapy. Recommend in-lab CPAP titration or APAP 5-15cmH2O.

Follow up with sleep specialist to ensure acclimation.



Howard Dean
23 Jun 2023

26.5
BMI

29.6
AHI

94
AVG SP0²

Good
Sleep^{co.}