

GoodSleep^{co}

Sleep Image Workflow




Getting started

→ Set up a test through Resdent

Use Resdent to set up a test for new or existing patient

 Resdent

 Patients

 Tests

 Lab Cases

 Referrals

 My Account

 Book a Call

 Logout

Patients

 Search patients

ALL

TEMPORARY THERAPY

PERMANENT THERAPY

Tests

Sync Recordings


+ Create Test

ALL

DATE ORDERED	SLEEP STUDY	PATIENT	NIGHTS ORDERED	STATUS
No tests found.				

Request Sleep Diagnostic

PATIENT



First Name

Last Name

Email

DOB

Phone

Gender

Create Patient

Select gender

→ Select Sleep Test type – Sleep Image sleep study or Sleep Image specialist reported

Sleep Image sleep study is auto-scored by Sleep Image and costs the clinician \$30 per night

A Specialist reported test provides a report and diagnoses from a Respiratory and Sleep Physician and costs clinician \$80 per night. This option will also allow for private health rebates should the patient require treatment for sleep disordered breathing

Request Sleep Diagnostic

PATIENT ☐ Create New

john smith - 01/03/23

x v

Sleep Study

Sleep Image Sleep Study

v

PATIENT ☐ Create New

john smith - 01/03/23

x v

Sleep Study

Sleep Image Specialist Reported

v

You must create a payment method or subscribe to ResDent to create a request

Add Payment Method

→ Select gender and number of nights then click 'Order Test'

Request Sleep Diagnostic

PATIENT ☐ Create New

john smith - 01/03/23

x v

Sleep Study

Sleep Image Sleep Study

v

Sleep Image

The SleepImage Rings stream data to the SleepImage Mobile App for upload, offering the highest level of comfort and convenience.

Want to learn more about Sleep Image? [Read More](#)

Nights

1 night

v

Gender

Male

v

You have 0 tests remaining before submitting this sleep study.

Order Sleep Test

GoodSleepco™

Sleep Image Workflow

→ The test will now show in your 'Tests' dashboard, including the patient key


The patient key will also be emailed to the patient. This is what will be used in the Sleep Image app by the patient to start the test.

DATE ORDERED	SLEEP STUDY	PATIENT	NIGHTS ORDERED
29 August 2023	Sleep Image	John Smith Sleep Image Patient Key: 5N7HNM	1

Set your patient up

Ring handover & App

- Sleep Image ring given to patient
- Patient to download Sleep Image app



→ Ensure patient has received email with the Sleep Image Patient Key

Clinician can also provide this from Resident

Patient - First night monitoring

→ Patient starts test

Patient puts the ring on and starts the test before bed using the Sleep Image app



→ Patient ends test when they wake up in the Sleep Image app

Keep the phone and ring close for a few minutes so the results can sync to the app

→ Email notification - Results ready

When the results are ready, you will be notified via email

→ Access the report in ResDent

You will be able to view these in the 'Tests' dashboard in Resdent (the specialist report may take a few days to appear but don't worry – you will get these back soon enough)

An interpretation guide for the test results will be supplied to you

					(sync recordings) / Create test
DATE ORDERED	SLEEP STUDY	PATIENT	NIGHTS ORDERED	STATUS	
4 August 2023	Sleep Image	[REDACTED]	1	Completed	<div>14 August 2023</div> <div>Endlink Test</div>



Sleep Image Workflow

Troubleshooting Tips

- **Ensure date of birth of patient is correct when creating the patient – this cannot be updated after**
- **Ensure there are no other Bluetooth devices in the room while doing the test**
- **Verify that you have a good internet connection by navigating to any webpage in your browser, for example:**
 - If you are on Wifi, try switching to cellular (or vice-versa)
 - If you have a VPN configured on your phone, that may cause issues. If you can disable it, try that.
 - If you have another mobile device in the household, such as a family members' phone or a tablet, try that and see if you are successful there.
- **If the recording time or quality is low once test results are received, please contact us hello@goodsleepco.health**



If you would like more information,
please contact us:

hello@goodsleepco.health