

capillus

BODY

User manual

Table of Contents

What is in the box	1
About the Capillus® BODY	2
The Capillus® BODY	2
Indications for use	2
Frequency of use.....	2
Duration of use.....	2
Warnings & Precautions.....	3,4,5
Safety Labels	5
Instructions for Use.....	6
Cleaning Instructions.....	7
Servicing.....	8
Storage Instructions.....	8
Disposal Instructions	9
Return Policy	9
Manufacturer's Guarantee	9
Risk of Loss, Other Terms of Sale.....	10
How to Contact Us.....	10
Device Specifications.....	11

What is in the box



1. **Capillus® BODY**

2. **Elastic Strap**

3. **Manual:**

A Physical Copy of the manual is available upon request at no additional cost. Call Customer service at (888) 571-5334

4. **USB-C Charging Cable with Power Adapter**

About the Capillus® BODY

The Capillus® BODY

Step into the science of light with the Capillus BODY belt, designed with dual-wave technology (630nm red light and 850nm near-infrared). Red light has been widely studied for its ability to stimulate mitochondrial activity, “the energy center of cells”, while near-infrared light penetrates more deeply, creating favorable conditions for circulation, tissue comfort, and recovery.

Engineered for Science. Designed for Comfort.

With powerful LEDs, a flexible design, and rechargeable convenience, the Capillus BODY makes advanced light technology simple to use at home or on the go. Whether after exercise, long workdays, or as part of your self-care, it's a modern approach to feeling restored, refreshed, and recharged.

Indications for use

The Capillus® BODY is an over the counter, at-home device. In just 15 minutes, the Capillus BODY delivers a gentle, warming experience that may help ease temporary aches, reduce everyday tension, and support the body's natural repair responses.

Research in the field of Low-Level Light Therapy (LLLT) has explored its role in promoting cellular energy, skin vitality, and favorable conditions for wound recovery—making it a powerful addition to modern wellness routines.

Frequency of use

For optimal results, use the Capillus® BODY 15 mins treatment once or twice a day. You may use the device on different parts of your body multiple times a day. To ensure the device is used optimally, allow the device to fully cool down before each use. Exceeding the recommended times or frequency may not increase the effect of the device.

Duration of use

The Capillus® BODY may be used indefinitely, provided it is operated according to its instructions.

Warnings & Precautions

Do not use the Capillus® BODY in:

- Conjunction with medications that are photo sensitizing.
- Products that are photosensitizing.

**Immediately stop using the device
if you experience any discomfort.**

Do not use the Capillus® BODY if you have any of the following medical conditions:

- Retinal and/or macular disease, ocular albinism, or any eye condition that causes sensitivity to light.
- Light sensitive skin disorder (e.g., lupus, porphyria).
- Light-induced headaches.
- Albinism.
- Melasma (i.e., condition that causes patches and spots).

Do not use the Capillus® BODY if you have any of the following conditions on the body.

- Suspicious lesions or skin cancer.
- Eczema.
- Psoriasis.
- Irregularly shaped, new, or changing brown spots.
- Broken skin, open wound, skin rash, or blisters.
- Tissue inflammation.
- Skin infection, cold sore, or warts.
- Active current outbreak of herpes simplex virus (cold sores).
- History of herpes outbreaks without receiving preventative treatment.
- Dermatological treatment where the skin has not fully healed, including neurotoxin, dermal filler, micro needling, laser, or chemical peel.

Warnings & Precautions Cont.

Use Warnings:

- Consult a medical professional before using if you are experiencing skin sensitivity or irritation.
- CAUTION: The light emitted from the Capillus® BODY may be harmful to the eyes. Do not stare at the light source.
- Do not use the Capillus® BODY if it has been left outdoors.
- To avoid the risk of temporary dry skin or hyperpigmentation, do not run back-to-back treatment sessions.

To reduce the risk of fire, thermal and chemical burns, electric shock, serious injury, or death, follow the precautions in this section and the procedures in this manual.

Charge Safely:





- Follow the manufacturer's "Instructions for use" to charge your Capillus® BODY.
- Only charge the Capillus® BODY with the provided USB-C charging cord and power adapter.
- Check if the voltage indicated on the equipment corresponds to the local mains voltage before you connect the equipment to avoid risk of electrocution or permanent damage to the equipment.
- Do not charge the Capillus® BODY in the bathroom or near water.
- Do not charge the Capillus® BODY after contact with water.
- Store the USB-C cord after charging, in a safe and dry place.

Warnings & Precautions Cont.

Supervision and Intended users:

- Do not use if any of the conditions in “Warnings & Precautions” (pages 3-4) apply to you.
- Do not leave the Capillus® BODY unattended when it is powered on.
- The Capillus® BODY is not intended for use by people with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instructions concerning use of the equipment. Please refer to the “Warnings & Precautions” before use.
- The Capillus® BODY is not intended for use by children under 18 years old. Children should be supervised to ensure they do not play with the Capillus® BODY.
- Small parts can be inhaled or swallowed by children—if any small parts break off or detach from the product, dispose of it safely.
- The Capillus® BODY is intended to be used by only one person.

Safety Labels

	Indicates the device should not be disposed of in the normal waste stream; it should be taken to a recycling center for proper disposal.
	Indicates Type BF applied part (60335); B=Body, F=Floating applied part. No conductive contact with USER.
	Indicates a medical device that needs to be protected from moisture.
	Indicates a battery that operates using lithium-ion.



Instructions for use

Please read the Warnings & Precautions and Instructions before first use.

Check the Product for Damage before each use.

- Inspect the Capillus® BODY for cleanliness, cracks, or sharp edges.
- Check the Capillus® BODY for signs of corrosion.
- These conditions can irritate or damage your skin.

1. Make sure your Capillus® BODY is fully charged. Please refer to the 'Charging Instructions' on how to charge your device.

2. To power on, press once and all LEDs will turn on. Press once again and all LEDs will turn off.

3. After you begin your session, the device is auto programmed to shut off after 15 minutes.

4. You may repeat the same procedure on different parts of the body. Allow the device to fully cool down before each use to ensure device optimal performance.

Charging Instructions

1. The CAPILLUS® BODY is equipped with an integrated rechargeable battery. To charge your CAPILLUS® BODY, connect it to the USB-C cord. Plug the USB-C cable into the power adapter, then plug the power adapter into a wall outlet.
2. Charging time: 1-4 hours.
3. During charging time, the LED indication will remain RED; it will turn off when charging is complete.

Cleaning Instructions

If it is necessary to clean the surfaces of your Capillus® BODY, please follow these instructions. Note: The Capillus® BODY has no parts that can be serviced by the user. If service is required, contact Capillus® at 1-888-571-5334.

MAKE SURE THE DEVICE IS TURNED OFF

UNPLUG THE CAPILLUS® BODY FROM ANY POWER CORD BEFORE CLEANING.

To clean the Capillus® BODY, only clean the fabric portion of the Capillus® BODY and not the LED directly.

Gently wipe the device with a damp cloth or use a soft brush to clean the elastic strap and the device. Allow it to air dry and make sure to keep it out of direct sunlight.

Do not use excessive pressure or force as this may damage the unit. Do not use acetone or any other solvents on any part of the device. Acetone and other solvents will damage the unit and void the warranty.

Keep It Dry.

- Keep the Capillus® BODY dry. If it gets wet, do not charge or use —the battery may be damaged. Follow disposal instructions.
- Do not use the BODY in wet surroundings (e.g., in a bathtub, shower, or swimming pool).
- Do not submerge the Capillus® BODY in water.
- Do not let water run into the Capillus® BODY.
- It is very important to keep the Capillus® BODY dry. Any contact with water could damage the device.
- Do not use or store the Capillus® BODY in the bathroom or near water.
- Do not charge or use the Capillus® BODY if it has been left outdoors.

Servicing

Protect the Battery from Damage.

- Do not try to open the Capillus® BODY. Do not replace the battery in the remote.
- Do not crush, drop, or damage the Capillus® BODY.
- Stop using or disconnect any power cable if the Capillus® BODY feels unusually hot or emits smoke or a bad smell.
- Follow "Disposal Instructions".

Do Not Repair or Modify.

- Do not attempt to repair or modify the Capillus® BODY.
- Do not crush or dispose of in fire.

Storage Instructions

Transport and Storage Conditions.

- Only transport and store the Capillus® BODY in the following environmental conditions: Keep between 0°C (32°F) and 25°C (77°F).
- Do not use the Capillus® BODY outdoors. Only use the Capillus® BODY indoors under the following environmental conditions Temperature: between 10°C and 30°C.
- A relative humidity range of 15% to 90% atmospheric pressure range of 700 Pa to 1060 Pa.

Disposal Instructions

Dispose of it in appropriate waste stream as per local ordinance and regulations.

Return Policy

All www.Capillus.com sales are final. No refunds or exchanges with the exception of the Manufacturer's Guarantee. If you wish to make a return as per above, contact Capillus® at 1 (786) 888-6249 to obtain a Return Material Authorization (RMA). Capillus® will not accept any shipments without an RMA. We recommend that you insure your shipment and register to receive tracking updates. Capillus® will only accept fully functional, like-new devices for returns. Capillus® will provide credit to your credit card or eligible refund via check.

If you received any products in a damaged condition, hang on to your shipment and DO NOT discard any damaged boxes, packaging, or products. Report the damage to us and we will determine if you should return the product for a replacement. You may be asked to provide pictures of your damage to help with the processing of your claim.

International orders will not receive any refunds for the shipping costs paid in the initial order. The refund will only be applicable to the price of the products.

Manufacturer's Guarantee

Capillus® warrants the device against defects in material or workmanship for a period of 60 days.

For more info visit: <https://www.Capillus.com/pages/terms-conditions>.

If you already have a Capillus® BODY. Do not forget to Register your Warranty at www.Capillus.com/warranty

If a defect in material or workmanship exists, Capillus® will repair or replace the necessary parts with new or reworked parts at no charge. Any replacement part will be covered under the original limited warranty. This Warranty is not transferable and does not cover damages caused by acts of God, misuse, negligence, accident, or modification of any part of the unit as determined by our service department. This Warranty does not cover damage due to improper operation or maintenance, connection to improper voltage supply, or attempted repair by anyone other than the company. Failure to use and maintain the device in accordance with the instructions in this manual will void the product warranty.

Please contact Capillus® customer service should you need to get warranty service or order replacement parts such as power adapter.

Do not fold, bend or disassemble any part of the Capillus® BODY and don't attempt to repair the device by yourself. Disassembling the device may damage it or may cause injury to you. If the Capillus® BODY is damaged, malfunctions, or comes in contact with liquid, contact Capillus®.

Risk of Loss; Other Terms of Sale

The risk of loss and title for items purchased by you pass to you upon our delivery of the items to the carrier. We do not take title to returned items until the item arrives at our fulfillment center. At our discretion, a refund may be issued without requiring a return. In this situation, we do not take title to the refunded item.

NO RESALE. Buyer shall use the device delivered to buyer for the buyer's user only, and shall not resell, transfer, exchange, or otherwise assign the device without first obtaining the prior written consent of Capillus®, which consent may be granted or withheld by Capillus® in its sole discretion.

How to Contact Us

You can contact us in writing at the following address:
Curallux, LLC 1715 NW 82 Ave Miami, FL 33126

By e-mail to our Capillus® Customer Team at: info@Capillus.com
By phone to our Capillus® Customer Team at: 1-888-571-5334



Device Specifications

Product Name:	Capillus® BODY
Model:	BODY-2025A
Rating Power:	60W
Device Dimensions:	400mm x 200mm x 4.25mm
Light Source Wavelength:	Red 630nm, Infrared 850nm
Timer:	15-mins sessions
Number of LEDs:	360 LED diodes (120 x 0.5W)
Weight:	0.58kg
Power Adapter:	5V/ 2.4A

Made in China

capillus

www.capillus.com

+1 (888) 571-5334