

osmo

MICROTONE USER MANUAL





NAME: MICRO CURRENT BEAUTY
INSTRUMENT MATERIAL: ABS PLASTIC
DIMENSIONS: 57X52X24 MM
WEIGHT: 35 GRAMS
BATTERY: 3.7V AT 170MA
FREQUENCY: 185 HZ
NOISE LEVEL: LESS THAN 50 DECIBELS
PATENT NUMBER: 202030103955.6

THE OSMO MICROTONE BRINGS PERSONAL TRAINING TO A NEW LEVEL OF BEAUTY. IT PROVIDES A SERIES OF CUSTOMIZED MICRO-CURRENT FACIAL EXERCISE MODES. THE MICRO-CURRENT CONTACT BALL CAN DIRECTLY CONDUCT MICRO-CURRENT TO THE SKIN TO HELP GENERATE COLLAGEN, REPAIR SKIN ELASTICITY AND TIGHTEN FACIAL CONTOURS.

HOW TO USE:

1. MAKE SURE YOU CLEANSE YOUR FACE THOROUGHLY REMOVING ANY EXCESS OIL. DO NOT USE AN OIL CLEANSER
2. APPLY SERUM ON YOUR FACE ESPECIALLY THE JAWLINE.
 - A. MAKE SURE BOTH THE DEVICE AND YOUR FACE / NECK ARE CLEAN AND DRY.
3. TURN ON THE DEVICE BY PRESSING DOWN THE POWER BUTTON FOR 2 SECONDS
4. CHOOSE THE INTENSITY BY PRESSING THE POWER BUTTON LIGHTLY. THE INTENSITY WILL BE SHOWN ON THE SIDE OF THE DEVICE.
 - A. EACH LIGHT MARKS THE INTENSITY OF THE MICRO CURRENT .
 - B. THE OSMO MICROTONE HAS THREE LEVELS
5. TONE & LIFT. MAKE SURE EVERY MOTION OF THE DEVICE IS OUTWARDS YOUR FACE IN AN UPWARD MOVEMENT

RECOMMENDED STROKES

1. 5 STROKES FROM THE BOTTOM OF THE JAWLINE TO YOUR EAR WITH YOUR JAWLINE INBETWEEN THE NODES
2. 5 STROKES FROM JUST UNDER YOUR NASAL FOLD TO YOUR EAR
3. 5 STROKES FROM THE LAUGHTERLINE UP TO YOUR TEMPLE, WITH YOUR CHEEKBONE INBETWEEN THE NODES
4. DO THIS ON EACH SIDE OF THE FACE APPLYING SERUM BEFORE TREATING EACH SIDE OF THE FACE.
5. FOR THE FOREHEAD, LIFT FROM YOUR BROW TO YOUR FOREHEAD IN AN UPWARD MOTION. COVERING THE WHOLE FOREHEAD
6. FOR THE NECK, STROKE FROM THE BOTTOM OF THE NECK TO YOUR JAWLINE. DO NOT USE IT IN THE THROAT AREA.

FEATURES:

1. RELIEVE MUSCLE TENSION: PRESS ACUPUNCTURE POINTS ON THE FACE TO RELIEVE MUSCLE TENSION.

2. MULTI-ACUPOINT STIMULATION: MULTI-ACUPOINT STIMULATION STIMULATES FACIAL REFLEX POINTS BY STIMULATING FACIAL ACUPOINTS TO CLEAR FACIAL MERIDIANS.

3. SCRAPER MICRO-CHAIN: THE SPECIALLY DESIGNED SCRAPER MICRO-CHAIN CAN FIT THE SKIN WELL AND HELP YOU MASSAGE YOUR FACE OR BODY TO ACHIEVE THE EFFECT OF FIRMING AND LIFTING.

4. USB CHARGING: LIGHTWEIGHT AND PORTABLE PACKAGING CAN BE PUT IN WALLET, HANDBAGS AND COSMETIC BAGS, SO IT CAN BE USED ANYTIME, ANYWHERE.

5. REDUCE SWELLING AROUND THE EYES: IT WILL REDUCE THE APPEARANCE OF SWELLING AND WRINKLES AROUND THE EYES, IMPROVE THE MICROCIRCULATION OF THE FACIAL BLOOD VESSELS, AND INCREASE THE FLOW OF BLOOD, LYMPH AND BODY FLUIDS. MAKING YOU LOOK PRESENTABLE, ENERGETIC AND HEALTHY.