

microsmooth STEP BY STEP



1. Thoroughly cleanse the treatment area with a drying cleanser. Rinse clean and pat dry with a towel.

2. Select the desired speed setting; low, medium or high speed by pressing the power button. Pressing the power button a 4th time will shut the device off.

3. Starting at the top of your face near the forehead hair line, begin making short feathery strokes as you work down and across the face and neck. You may gently glide over skin in



a treatment area 1-3 times positioning the microsmooth at a 45° angle. We **DO NOT** recommend using microsmooth over raised imperfections (including acne) the eyelid, directly on the nose or on the lips.

4. Apply a serum or moisturizer after treatment.

TIPS

1. **There's no wrong way to microsmooth** but aestheticians recommend gliding movements following the direction of use graphic for best possible results. Also, holding skin taut with opposite hand is recommended.




4. **The best time to microsmooth.** For best results, aestheticians recommend microdermabrasion in between dermaplaning treatments and when you aren't applying SPF or makeup directly after. This allows freshly exfoliated skin to receive maximum benefits of treatment. During your evening skincare routine is the ideal time to microsmooth.

5. **microsmooth and the sun.** Sonic dermaplaning and microdermabrasion is a rejuvenating exfoliating treatment. Take precaution with sun exposure before and after treatments.

2. **Areas to microsmooth include** the entire face and neck. Avoid any raised imperfections, the eye socket, nose and lips.

3. **The best way to hold microsmooth** is comfortably at a 45° angle.



6. **Sanitize your microsmooth after each use.** This ensures the best possible treatment results – we suggest using rubbing alcohol to clean the treatment tip.

7. **Frequency of treatments.** Dead skin and debris accumulate every 7-10 days. Spacing your Sonicsmooth/microsmooth treatments one week apart is ideal to maintain a smooth radiant complexion.

8. **Attaching and removing microsmooth treatment tips is easy.** To begin remove any older dermaplaning treatment tips, with the device powered off, grip the plastic harness around the edge and slide upward away from the device using your thumb. To attach a treatment tip slide the device handle through the opening along the tracks until microsmooth "clicks" into place.

NEED MORE INFO?

The entire operating manual including safety, care and warranty information is available at

www.sonicsmoothmanual.com

Contact customer service at 772-446-0149 or customerservice@michaeltoddbeauty.com

© 2022 MICHAEL TODD BEAUTY LP. All rights reserved.