


coldture

Red Light Therapy Face Mask

User Manual



Safety Instructions

For your safety and to ensure proper use, please read the user manual thoroughly before operating this device.

Caution:

Improper use may result in minor injury or damage to the device.

Notice:

Ensure all setup and operational procedures are followed as detailed in this manual to avoid potential harm or device damage.

Warning:

Avoid exposure to water or heat sources to prevent the risk of injury and device malfunction.


Specific Dangers

- Not suitable for users of pacemakers or similar devices.
- Keep away from water to prevent malfunction.
- Avoid use near heat sources like stoves or air conditioners.
- Do not use if damaged to avoid risk of electric shock.
- Keep out of reach of children.
- Do not use the mask for more than 20 minutes daily.

Setting Up Your Mask for First Use

1. **Prepare Your Skin:** Ensure your face is clean and dry.
2. **Fitting the Mask:** Place the mask on your face and adjust the strap for comfort.
3. **Turning On:** Press and hold the power button.
4. **Starting Treatment:** Press the mode button to select your desired therapy.
5. **Enjoy the Session:** Relax and adhere to the recommended therapy times.

How to Use the Mask




Turning On/Off

Press and hold the button for 2 seconds.

Recommended Timer Settings

Timer	Duration
T1	5 minutes
T2	10 minutes
T3	15 minutes

When first using the mask, use T1 for 5 minutes as a test.



Brightness Levels


For each mode, use Brightness on HIGH.

- Level 1: Low
- Level 2: Medium
- Level 3: High

Change Modes

Press the **M** button on the remote control to cycle through:

- **M1.** Glow Boost Therapy (Red + Near Infrared Light)
- **M2.** Clear Skin Therapy (Blue + Near Infrared Light)
- **M3.** Deep Rejuvenation Therapy (Red + Blue + Near Infrared Light)




Press the M (Mode) button

Remote Control On/Off

Press the **(I)** button to turn the remote control On and Off.

*Press the **(I)** button*



Light Therapy Modes and Customization

Adjust the timer to fit your skin's needs.

M1. Glow Boost Therapy (Red + Near Infrared Light)

- Benefits: Reduces wrinkles, fine lines, and boosts circulation.
- Recommended Use up to: 10 minutes, 4 times per week on high brightness.

M2. Clear Skin Therapy (Blue + Near Infrared Light)

- Benefits: Helps control acne, reduce inflammation, and redness.
- Recommended Use up to: 5-10 minutes, 4 times per week on high brightness.

M3. Deep Rejuvenation Therapy (Red + Blue + Near Infrared Light)

- Benefits: Deep skin penetration to reduce wrinkles, fine lines, and brighten skin.

Recommended Use up to:

- 10-minute sessions: Up to 4 times per week on high brightness.
- 15-minute sessions: Limit to 3 times per week on high brightness.


Note: For 15-minute sessions, start with three 10-minute sessions the previous week to condition your skin.

Adjustable Brightness Settings for Sensitive Skin

Our mask offers lower brightness settings for individuals whose skin does not take well to the higher power of the mask.

Charging the Mask

Connect to power using the USB-C cable.



Charging Indicators

- Red light: Charging
- Green light: Fully Charged

Additional Notes

- Charging time: 3 hours
- Battery life: Enjoy up to 40 minutes of use, allowing for 4 ten-minute sessions on low brightness or 3 ten-minute sessions on high brightness, per charge.
- Cannot be used while charging.

When the battery is low, the indicator light on the top of the mask will turn red. Please charge it in time to avoid affecting the next use of the product.

Charging in Controller

Please connect to power using the USB-C cable when the remote control is out of power.

Cleaning the Mask

Wipe with a damp microfiber cloth.

Do not use abrasive materials or strong cleaning agents.

Safety Certifications

Our device complies with CE, RoHS, FCC, GB, IE meeting international safety standards.

Optimal Results: Frequency and Duration of Use

For the best results LED Face Therapy Mask, follow these guidelines (do not use for more than 20 minutes per day):

M1. Glow Boost Therapy

- Recommended Use up to: 10 minutes at high brightness
- Frequency: 4 times per week

M2. Clear Skin Therapy

- Recommended Use up to: 5-10 minutes at high brightness
- Frequency: 4 times per week

M3. Deep Rejuvenation Therapy

Recommended Use:

- 10-minute sessions: Up to 4 times per week
- 15-minute sessions: Limit to 3 times per week

Note: For 15-minute sessions, start with three 10-minute sessions the previous week to condition your skin.

Voice Control

This product is voice-controlled, allowing you to interact with it and control the face mask through voice commands. The specific command terms and their corresponding action modes are as follows.

How to activate Mask: please say "Hi Lumi"

Voice Command	Mask Action
Turn on red light / Activate red mode / Start red therapy / Red light please / Begin red treatment	Start Mode 1, Red Light 630nm+850nm
Switch to red light / Change to red mode / Red light now / I want red light	Switch to Mode 1: Red Light 630nm + 850nm
Turn off red light / Stop red light / Red light off / End red treatment	Red light off and in standby mode
Turn on blue light / Activate blue mode / Start blue therapy / Blue light please	Start Mode 2, Red Light 460nm+850nm
Switch to blue light / Change to blue mode / Blue light now / I want blue light	Switch to Mode 2: Blue Light 460nm + 850nm
Turn off blue light / Stop blue light / Blue light off / End blue treatment	Blue light off and in standby mode
Turn on Magenta light / Activate Magenta mode / Start Magenta therapy / Magenta light please / Begin Magenta treatment	Start Mode 3, Red Light 460nm+630nm+850nm
Switch to Magenta light / Change to Magenta mode / Magenta light now / I want Magenta light	Switch to Mode 3: Magenta Light 460nm 630nm+ 850nm
Turn off Magenta light / Stop Magenta light / Magenta light off / End Magenta treatment	Magenta light off and in standby mode

Voice Control (continued)

Voice Command	Mask Action
Lowest brightness / Gentle mode / Dim light / Soft brightness / Level one brightness	Brightness: Low
Medium brightness / Normal mode / Standard brightness / Level two brightness / Balanced light	Brightness: Medium
Highest brightness / Strong mode / Maximum brightness / Level three brightness / Boost light	Brightness: High
Brighter please / Increase brightness / Level up brightness / More light / Step up brightness	Brightness +
Softer please / Decrease brightness / Level down brightness / Less light / Reduce brightness	Brightness -
Set 15 minutes / 15-minute session / 15 minutes please	Timer set to 15 minutes
Set 10 minutes / 10-minute session / 10 minutes please	Timer set to 10 minutes
Set 5 minutes / 5-minute session / 5 minutes please	Timer set to 5 minutes
Power off / Turn off / Shut down mask / End session	Face mask power off

Voice Control (Secondary Option)

Please note that voice control is provided as a secondary, backup option for this face mask. The performance of the voice recognition feature may be affected by various factors, including background noise and individual speech patterns.

For optimal performance:

- Please articulate your commands clearly and at a moderate pace.
- Ensure you are in a relatively quiet environment.

Important:

In the event that the device is unable to recognize your voice or respond to voice commands, please use the provided physical remote control as the primary method to operate the face mask.

Warranty & Guarantee Terms

Coldture offers a one-year guarantee covering material and craftsmanship defects under normal use, with exclusions for misuse, neglect, or unauthorized repairs.

Specification Details

Product Name	LED Face Therapy Mask
Product Size	640mm * 220mm * 5.8mm
Product Weight	375g
Light Wavelength	630nm, 460nm, 850nm
Output Power	6W
Chips Quantity	288pcs
Battery	1800mAh, use time 40 minutes
Function	Light Mode + Dimmable + Timer
Blue Light	460 ± 3nm
Red Light	665 ± 3nm
Near-Infrared Light	850 ± 3nm