

osmo

MICROLITE USER MANUAL





ELECTRIC NANOCRYSTALS BEAUTY INSTRUMENT

MICROLITE CHARGE AND PLUG IN ALL IN ONE MACHIN



1. Chip Installation
2. Dial to adjust the needle height
3. Gear Display
4. Long press the power supply for 1.5s to turn on the microcrystal switch
5. Power cord Android interface plug in instrument




The instrument should be waterproof, do not let liquid into the instrument. Cotton swabs are available after operation. Clean the inside of the instrument thimble to prevent the essence liquid from flowing into the instrument and causing jam.

GUIDE TO USE OF NANO ELECTRIC BEAUTY INSTRUMENT

APPLICATION INSTRUCTION GENERAL EDITION

- Step 1
Wash your face normally
- Step 2
Smear essence/other skin care products, Non alcoholic skin care product
- Step 3
Use the microneedle to make circles around your face, you can slide or tap



- Step 4
Apply a mask for 15-20 minutes using Non-alcohol no wash mask

- Step 5
Apply a cream to lock in nutrient moisture Non alcoholic no wash mask



Frequently Questions before use

hand and face should be thoroughly cleaned. Hands should be cleaned with hand lotion or soap. and face should be cleaned with facial cleanser or cleanser.

The chips are sterile and do not need to be sterilized.

Normal is 2 times a week, cuticle thin can be 1 times a week.

This recommended to import disposable independent packaging products with relatively single components, such as freeze-dried powder, raw liquid (hyaluronic acid,VC), water and light, moisturizing effect mask paste products.

Before introducing functional products, the product should be introduced in a small area behind the ear, and the skin will not have allergic symptoms such as redness, swelling, and itching Then apply it to your whole face.

Storage

1. Keep clean and dry, store in a well-ventilated environment, avoid high temperature exposure, damp, weight and impact.
2. Please do not put in the place where children can reach, away from children.
3. Please remove the nanocrystalline chip when storing the microcrystal.

Common problem in use

1. How to operate the import method?

When using, we should see the thickness and endurance of the cuticle of the skin, and the point of poor endurance, strong endurance .and cuticle Thick slips are generally used in areas such as the eyes, forehead, mouth triangle, cheek bones or neck because there is less subcutaneous fat make suggestions.

2. How to control timing adjustment and strength?

The time of first use is recommended to be less than 0.5, later can be adjusted according to their own situation, first gently touch the skin. Subject to skin feeling, No discomfort can be slightly increased strength.

3. How to take in various products?

Essence: after applying the essence on the face, apply it vertically to the skin with nanocrystals Face introduction, can also be smeared side import.

4. Can I use several products at once?

It is recommended to put no more than two products at a time, and the skin with poor tolerance can only import one product at a time.

Common problems after use

1. Do I need to apply other skin care products after importing?

It is recommended to use a moisturizing or repairing mask after importing. If the skin is dry, it can be combined with moisturizing cream or cream.

2. Can you touch the water? Can I wash my face and make up the next day?

It is recommended to import at night, do not touch the water on the day after import, if there is no discomfort symptoms such as redness, swelling, dry itching and so on the next day, it is ok Normal skin care and makeup, if uncomfortable symptoms, please use basic moisturizing products.

3. Can you enjoy the sun, spicy food, seafood and wine?

It is recommended not to be exposed to the sun within three days, eat less spicy, seafood, wine and other stimulating objects.

4. Can the chip be reused?

The chip is disposable and oxidizes when opened. After use, the chip sticks to the keratin and hair or the skin Bacteria breed and reuse is not recommended.

Usage steps

- 1, first clean the skin; Daily makeup removal, cleansing, skin toner can be used for secondary cleaning, remove the nano chip Mounted on the top of the micrometer.
2. Apply proper amount of original liquid or essence to your face.
3. Long press the power switch for 3 seconds to start the microcrystal and vertically touch the skin surface. It is recommended to stay at each contact point 3-5 seconds.
4. After the operation, please turn off the power and remove the nanocrystalline chip. Apply common skin care products. Cream, etc. (A wash- in mask before applying skin care products works best).
5. nanocrystalline chips can not be reused to avoid secondary pollution, please discard after use.

Operation method of microcrystal instrument

1. Pat and lift the face:

After the nanocrystalline wafer is installed, it is placed vertically on the skin surface, Slight indentation, speed, Dense: gentle up and down a lift, adjust the strength.

2. Sliding method:

After the nanocrystalline wafer is installed on the surface of the skin, slight depression can be visually observed in a straight line. While sliding, the nanocrystals should not leave the skin surface, and it is recommended to stay on the same part for 3-5 seconds No more than 2 times.

3. Circle method:

After the nano-chip is installed on the surface of the skin, the concave feeling can be visually detected, and the nano-chip can be absorbed in small circles at a constant speed The rice wafer should not leave the skin surface. It is recommended to stay on the same spot for 3-5 seconds and not more than 2 times.

Precautions

Notes for nano-electric microcrystal instrument:

1. Please do not remove or modify the charging device and microcrystal. It is strictly forbidden to use the charger when the power supply is damaged first. Otherwise, it can lead to electric shock. fire and charging damage.
2. Do not touch the USB charging cable with wet hands. If the charging cable touches any other liquid in fire or water. cut off the power immediately.
3. No short circuit charging device, no mechanical vibration or shaking instrument, no direct sunlight, no bathroom and other high temperature The product is not waterproof. it is strictly prohibited to rain or damp.
4. Please do not use the charge near the TV, radio and other electrical appliances.
5. Plug the head to focus on charging, pull the power cord will damage the charging line, charging must be microcrystal and charging line The use of uncertified charging cables may be hazardous in a ventilated area with a temperature of -10c C+55c C Breach of manufacturer's warranty.
6. Please keep the instrument out of the reach of children to avoid danger.

Skincare Precautions

1. Skin care products must be safe, effective, non-toxic and suitable for themselves.

Nanocrystalline Note:

1. If the disposable package is incomplete or broken, do not use it to avoid skin infection
2. Clean the operating parts and hands before operation.
3. For the first time, it is recommended to do a skin test on the inside of the wrist combined with skin care products, if itching occurs, etc.

Do not use when sensitive.

4. The strength of the chip should be moderate, and it is not necessary to press the chip hard or lift the wrist so that the chip can not touch evenly The skin.
5. About 10-20 minutes after using nanocrystalline chip, some people will have skin redness, which is a normal phenomenon If yes, the alarm will automatically disappear in 24 hours.
6. Except for the chip, other parts of the instrument are forbidden to contact with water and other liquids, so as not to damage the machine.
7. The same chip can only be used by one person at a time. It is strictly forbidden to share the same chip with others or disinfect the same chip After repeated use, otherwise it may cause skin allergy, red spots, inflammation and other adverse consequences.
8. When storing nanocrystals, discard the used wafers and store the unopened wafers and microcrystals in Keep well ventilated environment away from damp, pressure and impact, and keep away from children.
9. During use, the nanocrystalline wafer should keep moving and stay in the same position for 3-5 seconds, and the number of times should not exceed 2 times

Microneedle Precautions

1. Use of microneedles: Cleanse the face --- Apply the product to be introduced on the affected area --- Operate the microneedles for 10-15 minutes. It is recommended that when you perform the operation for the first time, you should try to be gentle and control the time to about 10 minutes (just stick it on the area to be done, no need to push hard). During the operation, the same part should not be used for more than 3 seconds, and the part to be done should be walked back and forth 3 times. Do not eat seafood within 3 days after microneedling.
2. Before using the product, it is recommended to apply a small amount on the skin behind the ears and use it only if there is no redness or itchiness.
3. After the operation, your face will be slightly red. Slight stinging and slight itching are normal phenomena. You can apply soothing ointment or ice mask to relieve pain.
4. The instrument should be waterproof and do not let liquid enter the instrument. Do not turn the instrument upside down during operation. Essence liquid can easily flow back from the needle into the tube body, making the machine unable to operate.
5. After using microneedle, do not touch tap water for 6 hours, as bacteria can easily cause infection. You can use saline, cool boiled water, or toner to wipe your face. Secondly, heavy makeup is not recommended.
6. After using microneedle, moisturize more, moisturize more, moisturize more, and protect yourself from the sun.