



Biella


Basic
Performance
Medical

Die Montage der Zubehörteile
führst du bitte wie folgt aus.






Rückfragen bitte an: support@wertach-quelle.com

Basic





- 1 GABRIEL-Tech Ferritkern
GFR9W88




Performance GABRIEL-TECH Edition



- 1 GABRIEL-Tech Ferritkern
GFR9W88
- 2 GABRIEL-Tech Chip
GT-WQ-DEGDW90BL44 - dunkelblau
- 3 GABRIEL-Tech Chip
GT-SI-DEGDL90SI48 - hellblau/silber
- 4 GABRIEL-Tech Chip
GT-WQ-DEGDM60SC88 - schwarz



Medical UMH Edition



- 1 GABRIEL-Tech Ferritkern
GFR9W88
- 2 GABRIEL-Tech Chip
GT-WQ-DEGDW90BL44 - dunkelblau
- 3 GABRIEL-Tech Chip
GT-SI-DEGDL90SI48 - hellblau/silber
- 4 GABRIEL-Tech Chip
GT-WQ-DEGDM60SC88 - schwarz
- 5 UMH GOLD/RHODIUM light
Hinweis: Das UMH Modul erst nach
der Erstsplüfung der Filter einsetzen.

